

INNTEG Hydrogen Generators

High Purity & Ultra High Purity

Need a **continuous supply** of **high purity Hydrogen**?

Need a **safe and reliable** Hydrogen source?

Need a Hydrogen Generator for your **gas chromatography**?

INNTEG Hydrogen Generators are your **PERFECT** choice



- Stylish Design
- Modern System Technology
- No Maintenance
- Best Price Performance

INNTEG

Hydrogen Generators

Specification		
Model	HP-HG250	UHP-HG250
Technology	Proton Exchange Membrane (PEM)	
Purification	Perma Pure Membrane Drying Tubes	Perma Pure Membrane Drying Tubes + Heatless Regenerative Desiccant Dryer
Flow Rate	250 mL/min	
Purity	>99.9996%	>99.99996%
Pressure	0.5 to 11 barg	1.4 to 11 barg
Safety	Internal and External Leak Detectors	
Recommended Application	Fuel Gas (FID, NPD & FPD)	Fuel Gas (FID, NPD & FPD) and Carrier Gas

Please contact us for more details.

 (852) 3910 1533

 hkcsr@tegent.com.hk